

**MEDICAL MICROBIOLOGY BRANCH
VIROLOGY SECTION
SPECIMEN REQUIREMENTS FOR SERUM SPECIMENS SUBMITTED FOR
DENGUE IGM CAPTURE ENZYME LINKED IMMUNOSORBENT ASSAY
(ELISA)**

The controlled version of this document is the electronic version located in the folders described in section 5.0 of the Medical Microbiology Branch Document Control Standard Operating Procedure. If this is a printed copy, it is an uncontrolled version and must be verified with the controlled version BEFORE use.

SPECIMEN REQUIREMENTS FOR SERUM SPECIMENS SUBMITTED FOR DENGUE IGM CAPTURE ENZYME LINKED IMMUNOSORBENT ASSAY (ELISA)

Methodology:	IgM Capture Enzyme Linked Immunosorbent Assay (ELISA)
Performed:	<p>The ELISA for IgM antibodies is the most useful serologic procedure available for determining recent infections by dengue virus.</p> <p>Serological results must always be interpreted in the context of compatible signs and epidemiologic data.</p>
Criteria for testing:	<p>Specimens meeting the case definition set by the Disease Investigation Branch (DIB) of the Disease Outbreak Control Division (DOCD) of the Department of Health. Approval by the DIB/DOCD, is <u>required</u> before specimens are tested.</p> <p>Clinical signs and symptoms compatible with Dengue fever.</p>
Turn-Around-Time:	<ul style="list-style-type: none">• Results are reported out 48 -72 hours after testing.• Priority testing<ol style="list-style-type: none">1) May be requested by the Disease Investigation Branch on a case by case basis.2) Prior notification of the laboratory is required before submission of the specimens.
Specimen (type) required:	A minimum of one (1) mL of serum.
Specimen Collection:	Serum collected as early as possible after onset of illness and placed in a sterile container. Most patients become positive by the sixth day after onset of illness. Multiple specimens may need to be collected in case the specimen is collected too early. Detectable IgM antibody usually persists for a maximum of 60 to 90 days.

Specimen storage, packing and transport:	Ship specimens in sterile, screw-capped tubes. Follow instructions for Class B – Biological Substance of the U.S. Department of Transportation (U.S. DOT) and International Air Transport Association (IATA) for packing and shipping.
Specimen submission:	Submitters (Clinical Laboratories and Epidemiology Specialists of the Disease Outbreak Control Division) should notify DIB and the laboratory prior to the submission of specimens.
Requisition Form:	Submitter is responsible for completing MMB Form 81.3 (including but not limited to the following information: patient identifier, date of onset of illness, clinical signs and symptoms, travel history, immunization history, name and address of submitter).
Unacceptable conditions:	<ul style="list-style-type: none"> • Specimen that is leaking; • Improper container or handling; • Obvious microbial contamination; • Specimen quantity is insufficient to perform the test; • Unlabeled specimens; • Improper or incomplete requisition form (MMB Form 81.3); • No requisition form; • Specimen label does not match the requisition.
Stability:	Serum specimens may be stored at 2°C to 8°C for a maximum of 48 hours. For longer storage, freeze at ≤ -20°C.
Normal Value:	No IgM antibody detected to Dengue Virus.
Limitations:	<p>Most serologic tests cross react to a variable extent with other members of the flavivirus family. This is a particular problem in areas where dengue and Japanese encephalitis or yellow fever occurs together.</p> <p>Dengue virus reaches peak titers in the blood before acute illness. Most patients become positive by the sixth day after onset of illness. Detectable IgM</p>

antibody usually persists for a maximum of 60 to 90 days, and is a good measure of recent infection.

Result Notification: Laboratory results are reported to the submitters and Disease Investigation Branch of the DOH Disease Outbreak Control Division.

Test performed at: Virology Section, Medical Microbiology Branch
State Laboratories Division
Department of Health
2725 Waimano Home Road, 2nd Floor
Pearl City, Hawaii 96782

Contact: Roland Lee, Supervisor, Virology Section
(808) 453-6705

Lynne Washio, Supervisor, Viral Isolation and
Detection Unit
(808) 453-6705

Dayna Ornellas, Chief, Medical Microbiology Branch
(808) 453-6700

Validated By:



Roland H.K.C. Lee
Supervisor, Virology Section

10/22/25

Date

Reviewed By:



Dayna Ornellas, M.S.
Chief, Medical Microbiology Branch

10/22/25

Date

Approved By:



Edward P. Desmond, Ph.D.
Administrator, State Laboratories Division

OCT 21 2025

Date